




Sartorius Enclosed Scale  
Model: CPA324S




**7**  
The door to the glass enclosure flips back on the left side. Please make sure it's closed to avoid inaccuracies.




**8**  
Unit power cord is located on the rear.



**9**  
When not in use, please recover using plastic slip cover.



**11**  
Leveling bubble



**10**  
Adjustable feet for leveling surface.



**6** Tare button  
**1** Unit power button  
**2** Cancel/Clear button  
**3** Calibration button  
**4** Function button  
**5** Data Output button

# SETUP

A

Ensure unit is plugged in and area is free of debris. **8**

B

By rotating feet, adjust leveling bubble so that bubble is in center. **10**

C

Once leveled, ensure door is closed and turn unit on. **1**

D

Wait while the unit performs a display test.

E

Press the Tare button to zero out the balance. **6**

F

Place sample on the balance.

# OPERATION

A

With setup complete, place samples and record data.

B

To erase a weight/start over, you may press cancel or tare. **2**

C

For advanced functions such as calibration, see user's manual. **3**

D

To print out data, press the data output button. **5**

E

When finished, press the power button to turn off. **1**

F

Please replace plastic slip cover when done. **9**

# HAZARDS



**ELECTRICAL HAZARD**  
This equipment contains electrical units that may come into contact with skin or liquids.



**SLIP HAZARD**  
Because chemicals are being moved near this unit, a slip hazard may be created.



**PINCH HAZARD**  
The glass door to the balance may shut and cause a pinch hazard.



**CHEMICAL HAZARD**  
Because chemicals are being moved to the balance surface, chemical spills are possible.



**FIRE HAZARD**  
Any improper use of this unit may cause a fire hazard.

# PREVENTION



**EYE PROTECTION**  
Wear safety glasses while operating.



**PROTECTIVE GLOVES**  
Wear appropriate gloves for heat and chemical usage.



**NON-SLIP, CLOSED SHOES**  
Wear non-slip, closed shoes to avoid spills.



**READ MANUAL**  
Become familiar with this guide before operating.

For additional help contact the Lab Space Manager