


Time-Resolved, *In-Situ* Neutron Diffraction Studies of Hydrothermal Reactions



ACNS, Knoxville, TN - June 24th 2002

Dermot O'Hare,
Inorganic Chemistry Laboratory
University of Oxford, UK

Solid State Reactions that we have Investigated using Time-Resolved *In Situ* Powder Diffraction


Synchrotron X-Ray Experiments

- Crystallisation from Solution
- Intercalation Reactions
- Hydrothermal Reactions
- Sol-Gel Syntheses
- Solid State Syntheses in Molten-Salt Fluxes
- Classical Solid Phase Reactions

Neutron Experiments

- Hydrothermal Reactions

Hydrothermal Synthesis



Many experimental variables:

- starting mixture composition
- choice of reagents
- temperature
- time
- % fill of reaction vessel

Sealed reaction vessel

- observation of reaction course is difficult

It remains very difficult (impossible) to predict the outcome of new reactions.

⇒ need to develop methods for *in-situ* study

Key Features of Energy Dispersive X-ray Powder Diffraction using Synchrotron Radiation

Advantages

- Very fast data collection times; typically (1 – 30 sec)
- High energy incident X-ray photon beam; give excellent penetration through cell materials
- Large volume cells; allows us to reproduce laboratory scale syntheses


Disadvantages

- Incident beam profile not well characterised
- Poor resolution, especially for high d-spacings
- Structure solution and/or refinement impossible
- Strongly absorbing samples

Time-Resolved *In-Situ* Powder Neutron Diffraction

- Cell Design
- Gem Diffractometer
- Hydrothermal Chemistry
 - Synthesis of Zeolites
 - Synthesis of t -BaTiO₃


Schematic Diagram of the Oxford/Isis Hydrothermal Cell



- Constructed from null-scattering Ti-Zr alloy
- 25 mL volume
- 4 cm high window of material exposed to beam
- Internal surfaces sputtered with gold


Walton *et al.* *Rev. Sci. Instrum.* **70 (1999) 3391**

Photograph of the Oxford/ISIS Hydrothermal Cell




Hydrothermal Reaction Cell

- constructed from null-scattering Ti-Zr alloy
- 25 mL volume
- 4 cm high window of material exposed to beam
- Internal walls sputtered with 100 μ m of gold



Close-up Photograph of the Assembled Hydrothermal Cell

- Null Scattering Alloy
- Copper Heating Blocks require BN Shielding
- Facility for stirring reactions



Photograph of the Fully Assembled Hydrothermal Cell




GEM: Shining Bright



Instrument Parameters

DETECTORS:	4200 \uparrow 6600 “Conventional” 7-bank structure (resolution-focused).
SOLID ANGLE:	2.85 \uparrow 4.2 Sterad.
CHOPPERS:	3 choppers (λ -selector, frame overlap, nimonic).
Q-RANGE:	Very high Q_{MAX} ($E \leq 1$ eV, $l \geq 0.28$ Å).
RESOLUTION:	~0.2% in back-scattering.
DATA SET SIZE:	64 Mbytes \uparrow 128 Mbytes
DATA RATE (PEAK):	6 Gbytes-500 files/day
DATA RATE (AVE):	1 Gbyte-40 files/day


Schematic Diagram of the Detector Array on the New GEM Diffractometer at the Isis Facility





22/2/2000


Hydrothermal Cell Background



The cell gives rise to virtually no background and provides a chemically inert environment.

Rietveld Refinement of the TOF Neutron Powder Diffraction Data for BaTiO₃

In the Hydrothermal Cell at 120 °C



Refined Parameters

- Background, peak profile
- Unit cell parameters of BaTiO₃ and TiO₂

Time-Resolved *In-Situ* Powder Neutron Diffraction


- **Hydrothermal Chemistry**

Synthesis of Zeolites

Synthesis of *t*-BaTiO₃


Hydrothermal Synthesis of Sodalite

Synthesis:




- Apertures formed by 6-rings only
- Cubic, $\overline{\text{P}}43\text{n}$
- $a = 8.9 \text{ \AA}$

Crystallisation of Zeolite A and Hydroxosodalite




Low NaOH
concentration



Sodium zeolite A (LTA)



High NaOH
concentration




Hydroxosodalite (SOD)




Long reaction
times

Time-Resolved Powder Neutron Diffraction Data for the Hydrothermal Transformation of Zeolite A to Hydroxosodalite


at higher pH:
collapses into the denser sodalite phase



direct crystallisation of sodalite




140 °C:





Growth Kinetics Obtained from *in situ* Powder Neutron Diffraction Data

$\text{Al}_2\text{O}_3 : 2\text{SiO}_2 :$
 $3.5 \text{ NaOD} : 20 \text{ D}_2\text{O}$



- Follows Ostwald's rule of successive crystallisations

Hydrothermal Synthesis of BaTiO₃



t-BaTiO₃

Hydrothermal method provides;

- slow temperature synthesis
- small (submicron) particle size
- good particle size distribution

- **Ferroelectric material**
- **High permittivity**
- **Applications**
DRAM, electro-optic devices, capacitors




Following the Hydrothermal Crystallisation of BaTiO_3 using *in-situ* Powder Neutron Diffraction

POLARIS Diffractometer

12 hours: BaTiO_3

125 °C



15 minutes per spectrum



15 minutes: $\text{Ba}(\text{OD})_2 \cdot 8\text{D}_2\text{O} + \text{TiO}_2$

Walton *et al.* Chem. Commun. (2000) 1267

Growth Kinetics of BaTiO₃ from *in-situ* Powder Neutron Diffraction




First *in situ* study of the hydrothermal crystallisation of BaTiO₃

- Large amount of TiO₂ dissolves before BaTiO₃ is observed
- Small particle-size amorphous Ti source - markedly faster crystallisation

⇒ dissolution-precipitation mechanism

- Suggests a reaction between solution phase Ti(OH)_x^{4-x} and Ba²⁺ ions


Crystallisation at Barium Titanate



Each of the tof neutron data sets was measured in 5 minutes on GEM


Time-Resolved in-situ Powder Neutron Diffraction Data for the Hydrothermal Crystallisation of BaTiO₃

Synthesis Conditions: 1.1BaCl₂: TiO₂: 3.33 NaOD: 16.7 D₂O at 200 °C



Normalised Integrated Bragg Peak

Synthesis Conditions: 1.1BaCl₂: TiO₂: 3.33 NaOD: 16.7 D₂O at 200 °C




Crystallisation of Barium Titanate




125 °C

- Competitive formation of a second phase
- Identified as Ba_2TiO_4




Time-Resolved in-situ Powder Neutron Diffraction Data for the Hydrothermal Crystallisation of BaTiO₃

Synthesis Conditions: 1.1BaCl₂: TiO₂: 3.33 NaOD: 16.7 D₂O at 125 °C




Rate of Crystallisation of BaTiO₃

Effect of temperature



Rate of Crystallisation of BaTiO_3

Effect of NaOD Concentration



Nucleation and growth is
independent of pH

Conclusions from this Study of the Hydrothermal Synthesis of *t*-BaTiO₃

- BaTiO₃ begins to crystallise at *ca.* 20 mins and is complete at *ca.* 40 mins.
- Both Ba(OH)₂ and BaCl₂ starting materials dissolve before the onset of crystallisation.
- A significant amount of TiO₂ dissolves immediately on heating, and before the onset of BaTiO₃ crystallisation.
- At 125 °C a transient crystalline intermediate (Ba₂TiO₄) is observed.
- The nucleation time and growth rate is independent of pH.

Conclusions

In situ diffraction methods provide an efficient means of following hydrothermal reactions under *real conditions*

Energy-Dispersive X-ray Powder Diffraction

- high time resolution (< 1 minute)
- kinetic data of unrivalled quality
- modelling kinetic data gives mechanistic information
- transient crystalline phases can be observed

Powder Neutron Diffraction

- diffraction data have high resolution
- structural information from *in situ* diffraction data
- higher neutron fluxes and the use of position-sensitive detector banks will vastly improve data acquisition rates

Acknowledgements

Time-Resolved Experiments

Dr R. Walton
Dr F. Millange
Dr A. Norquist

Collaborators

Prof G. Férey (Versaille)
Dr T. Loiseau (Versaille)
Prof C.N.R. Rao (Bangalore)

Instrument Scientists

Dr R. Smith (ISIS)
Dr P. Radaelli (ISIS)
Dr T. Hanson (ILL)

Financial Support

EPSRC
Royal Society
CLRC Laboratories