

*“Any sufficiently advanced technology is indistinguishable
from magic”*

Arthur C. Clarke


New Tricks with an Old Tool:
Neutron Spin Echo
Applied to
Reflectometry and SANS

By


Roger Pynn

Los Alamos National Laboratory
and the University of California at Santa Barbara


Nanoscience & Biology Need Structural Probes for 1-1000 nm




10 nm holes in PMMA




CdSe nanoparticles




Peptide-amphiphile nanofiber




Actin



Si colloidal crystal



Structures over many length scales in self-assembly of ZnS and cloned viruses



Thin copolymer films

Current Limitations of Neutron Scattering for Nanoscience and Biology


- Difficult to probe lengths much larger than a 100 nm with SANS
- Measurement of reflectivity is often limited by diffuse scattering or incoherent background
- No routine measurement of in-plane structures of thin films
- Time scales for kinetic measurements with SANS are quite long so that special systems have to be chosen for such measurements
- Very limited measurement of inelastic scattering from thin films

Tight collimation is incompatible with high signal intensity

Neutron Spin Echo (NSE) uses Larmor Precession to “Code” Neutron Velocities


- A neutron spin precesses at the Larmor frequency in a magnetic field, B . $\omega_L = \gamma B$
- The total precession angle of the spin, ϕ , depends on the time the neutron spends in the field

$$\phi = \omega_L t = \gamma B d / v$$



$$\text{Number of turns} = \frac{1}{135.65} \cdot B[\text{Gauss}] \cdot d[\text{cm}] \cdot \lambda[\text{Angstroms}]$$

In NSE, Neutron Spins Precess Before and After Scattering & a Polarization Echo is Obtained if Scattering is Elastic




By “Inclining” the Magnetic-Field Region, Spin Precession Can Be Used to Code a Specific Component of the Neutron Wavevector

If a neutron passes through a rectangular field region at an angle, its total precession phase will depend only on k_{\perp} .

$$\omega_L = \gamma B$$


$$\phi = \omega_L t = \gamma B \frac{d}{v \sin \chi} = \frac{KBd}{k_{\perp}}$$

with $K = 0.291 \text{ (Gauss.cm.}\text{\AA)}^{-1}$



With $Bd = 1000 \text{ Gauss.cm}$ & $k = 1.5 \text{ \AA}^{-1}$ we get 1 radian change in ϕ for $\delta\chi \sim 0.2^\circ$ at $\chi = 45^\circ$ or for $\delta\chi \sim 0.01^\circ$ at $\chi = 10^\circ$

A Simple Example of Tilted Fields: SANS



- Any unscattered neutron ($\theta=0$) experiences the same precession angles (ϕ_1 and ϕ_2) before and after scattering, *whatever its angle of incidence*
- Precession angles are different for scattered neutrons

$$\phi_1 = \frac{KBd}{k \sin \chi} \text{ and } \phi_2 = \frac{KBd}{k \sin(\chi + \theta)} \Rightarrow \cos(\phi_1 - \phi_2) \approx \cos \left[\frac{KBd \cos \chi}{k \sin^2 \chi} \theta \right]$$

$$P = \int dQ . S(Q) . \cos \left[\frac{KBd \cos \chi}{k^2 \sin^2 \chi} Q \right]$$

Polarization proportional to
Fourier Transform of $S(Q)$


Spin Echo Length, $r = KBd \cos \chi / (k \sin \chi)^2$

How Large is the Spin Echo Length for SANS?


$Bd/\sin\chi$ (Gauss.cm)	λ (Angstroms)	χ (degrees)	r (Angstroms)
3,000	4	20	1,000
5,000	4	20	1,500
5,000	6	20	3,500
5,000	6	10	7,500

- It is relatively straightforward to probe length scales of ~ 1 micron
- Note amplification of r by decreasing χ and increasing λ


Discriminating Between Specular and Diffuse Neutron Reflection with Better than 0.5 mrad Angular Resolution






Rotator and precession films
were $\sim 30 \mu\text{m}$ thick Permalloy
($\text{Ni}_{0.8}\text{Fe}_{0.2}$) easily saturated in 10 Oe



Experiment at Hahn-Meitner-Institut





A 30 μ Permalloy Film used as a Neutron Spin Rotator




Experimental Results for Reflectometry

- 30 μm permalloy films magnetized by a guide field work well as $\pi/2$ rotators
- 30 μm permalloy films, magnetized by a guide field, work as “precession fields” but thickness variation limits echo amplitude to ~ 0.5
- Change in the echo amplitude was observed when beam was diffusely scattered
 - Our apparatus was sensitive to changes of about 0.02° in neutron trajectory




Using “Sign Reversal” to Implement Angle Coding

An element that performs a π rotation about an axis in the precession plane changes the sign of prior precession angles



- Total net precession angle = $(\phi_4 + \phi_5 + \phi_6) - (\phi_1 + \phi_2 + \phi_3) + 2(\phi_2 - \phi_5)$
- The first two terms depend on neutron velocity only, last term depends on velocity *and* angle of neutron trajectory
 - Requires suitably oriented planar π rotators (flippers)
- Can be set up to encode two, mutually perpendicular, trajectory angles

Using Sign-Reversal for Reflectometry



The NSE Technique May Provide a Way to Enhance Signal Intensity & Resolution for *Elastic* Scattering

- The method can be applied to achieve good resolution without the need for tight monochromatization and collimation
 - Extend size range for SANS
 - ~100x gain in measurement speed for SANS at same resolution & accuracy
 - Separate specular and diffuse (or incoherent) scattering in reflectometry
 - Measure in-plane ordering in thin films (Felcher)
- The method can be combined with real-space focusing techniques such as focusing supermirrors or lenses
- Thin permalloy films work as neutron spin rotation devices and allow implementation of “sign reversal” method
 - Only demonstrated for constant wavelength spectrometers but can probably be achieved for white beams