

**How GDA helps to turn raw data into results
Or at least lets you decide
if you are collecting reasonable data**

Tobias Richter
for the Diamond data acquisition & scientific software teams



NOBUGS 2010 Gatlinburg



Introduction



@ Diamond

- is the data acquisition system used on (almost) all beamlines
- has integrated visualisation and analysis features
- caters for a wide range of disciplines
- written in Java
- has a Jython scripting interface
- is a client/server based system
- GUI client written using the Eclipse RCP framework
- is open source software, see www.opengda.org

Integrated Data Acquisition and Processing

- give the user control over the data he is collecting by displaying in a meaningful way
- do not waste beamtime on bad samples
- user is more likely analyse the data in more detail and publish
- can be used for decision making in automated experiments
- should not interfere with data collection

Integrated Data Acquisition and Processing


- give the user control over the data he is collecting by displaying in a meaningful way
- do not waste beamtime on bad samples
- user is more likely analyse the data in more detail and publish
- can be used for decision making in automated experiments
- should not interfere with data collection
- collection and processing at the same time?

Integrated Data Acquisition and Processing



Exafs Analysis

- GUI based
- no processing data is written to the file
- scans with scalar data ($h\nu$ vs. I_0 and I_t) – low data rate
- background subtraction, first derivative and fourier transformation are performed



Exafs setup

Data Acquisition Client - Beamline BL188

File Edit Window Help

Plot Single Scan Multiple Scan

XAS_Parameters.xml B18Sample_Parameters.xml Detector_Parameters.xml Output_Parameters.xml

XAS Parameters

Element Pt
Edge L3
Edge Energy 11564 eV
Core Hole 5.31 eV

Scan Parameters

Initial Energy 11304 eV
Final Energy 13250 eV 21.036 Å⁻¹
Edge Region Gaf1/Gaf2
Gaf1 25
Gaf2 10
A 11404.7 eV
B 11510.9 eV

Step Parameters

Pre-Edge Energy Step 5 eV
Pre-Edge Time Step 0.5 s
Edge Energy Step 1 eV
Edge Time Step 0.5 s
Exafs Step Type E
Exafs Step Energy 2 eV
Exafs Time Type Constant Time
Exafs Time Step 0.5 s

Number of scan points 981 points
Estimated time 00:22:50
EXAFS Graph


XAS Scan XML

Single Scan Controls


Pause Run Stop View Scan

zj121856 Baton held N...g 386M of 565M Not logged in

Exafs data



Exafs fourier




Small Angle Scattering

- processing deployed on I22
which serves physicists, chemists and biologists
- experiments collect from 1MB to 40GB
and run in 200ms to 4h
- processing runs on server
- processing parameters and results are written to file along
with raw data

Automatic Reduction

- normalise
- detector response
- background subtraction
- average
- invariant
- sector integration



GDA standard tools to visualise and investigate data are used to define reduction parameters (i.e. sector integration)


Data Processing Features

- delivers I vs. q
- reduction results saved to NeXus with raw data
- framework allows any data operation
- extensible (also from scripting) and configurable on the fly
- elements can be combined in a chain or tree
- all data (reduced and raw) for (at least) one frameset needs to fit into memory

Implementation Details

- detectors present themselves as a detector system with a configureable set of virtual detectors
- each detector (physical or virtual) contributes a partial NeXus tree to the savable data
- each detector receives a handle to the tree structure to read and add data
- when the physical detectors are read, the scan continues
- processing detectors can keep state, if they want to operate across scan point boundaries

Threading



The number of processing threads can be configured at runtime.

Problem

Users need to give up the first few minutes of their beamtime for measuring background, calibrants, etc.

Ongoing Work

- configurable visualisation in domain specific plots (Kratky, Guinier, Porod, Zimm, etc)
- re-evaluation for offline analysis
- simple analysis: peak fitting / finding / following

Outlook

- tuning for performance responsiveness for more exotic cases
- more advanced analysis: domain specific configurable workflow run outside of GDA (using EDNA), results displayed in GDA
- decision making
- deployment on other beamlines (at Diamond and elsewhere)

A collage featuring a photograph of a combine harvester in a field and abstract colorful flowing lines. The photograph shows a combine harvester with a red wheel and a blue body, harvesting golden-brown grain. The abstract lines are multi-colored (yellow, green, blue, pink) and flow across the white background, partially overlapping the photograph.

