



ORNL




NEUTRON SCIENCES




**Eppendorf MiniSpin Plus
Micro Centrifuge**



SETUP



OPERATION



HAZARDS



ELECTRICAL HAZARD

This equipment contains liquid near electrical connections.



PROJECTILE HAZARD

Use of an imbalanced or damaged rotor may result in a projectile hazard.



PINCH HAZARD

This equipment contains moving parts near the lid/agitator that may cause physical harm.



INHALATION HAZARD

Improperly sealed tubes or spilled liquids may atomize or vaporize in the chamber.

PREVENTION



EYE PROTECTION

Wear safety glasses while operating.



PROTECTIVE GLOVES

Wear appropriate gloves for heat and chemical usage.



NON-SLIP, CLOSED SHOES

Wear non-slip, closed shoes to avoid spills.



READ MANUAL

Become familiar with this guide before operating.

For additional help contact the Lab Space Manager