



**ORNL**




NEUTRON SCIENCES




**Eppendorf 5415D/5415R  
Microcentrifuge**



# SETUP



# OPERATION



# HAZARDS



## ELECTRICAL HAZARD

This equipment contains liquid near electrical connections.



## PROJECTILE HAZARD

Use of an imbalanced or damaged rotor may result in a projectile hazard.



## PINCH HAZARD

This equipment contains moving parts near the lid/agitator that may cause physical harm.



## INHALATION HAZARD

Improperly sealed tubes or spilled liquids may atomize or vaporize in the chamber.

# PREVENTION



## EYE PROTECTION

Wear safety glasses while operating.



## PROTECTIVE GLOVES

Wear appropriate gloves for heat and chemical usage.



## NON-SLIP, CLOSED SHOES

Wear non-slip, closed shoes to avoid spills.



## READ MANUAL

Become familiar with this guide before operating.

For additional help contact the Lab Space Manager