

Science Case for Better Neutron Detectors for Scattering

Lowell Crow
Spallation Neutron Source
Experimental Facilities Division

May 29-30, 2003

Spallation Neutron Source at Oak Ridge



SNS – Complete in 2006



Why Neutrons?



Neutron Interactions: coherent scattering
incoherent scattering
magnetic scattering
capture

Thermal Neutrons:

Wavelengths comparable to atomic spacing

Velocities comparable to atomic vibrations

Scattering cross section is nuclear, not electronic

Distinguish isotopes, hydrogen

Sensitive to magnetic moments

Neutrons are a bulk probe – interactions are weak

Detector performance parameters



Neutron capture efficiency at maximum operating energy

50% at 50 meV to 50 eV

Counting efficiency

near 100%

Position resolution

2 cm to 100 μm

Gamma-ray rejection

10^{-6} to 10^{-8}

Rate

10^2 to 10^6 n/s/pixel

10^5 to 10^8 n/s (whole area)

Time resolution (at a known flight path distance)

1 μs to 10 μs

Diffraction

Elastic scattering yields a diffraction pattern dependent on the atomic arrangement of the sample.


Instruments are designed for the type of sample, or for the specialized environment.

Powder

Single Crystal

Liquid/Amorphous/Short-range ordered

Methane
Hydrate



Loveday et al.

High pressure

High magnetic field


Engineering materials

Residual Stress


Texture

Some Diffraction Instruments


Wollan and Shull at X-10



GEM
at
ISIS



Vulcan at SNS



Powder Diffraction



Nuclear scattering gives atom positions, especially in materials with a combination of heavy and light atoms

- Metal oxides – *oxygen positions*
- High temperature superconductors
- Parametric studies
 - Temperature dependence
 - Processing dependence

Magnetic scattering -- main source of information about moment ordering in magnets


1st observation of antiferromagnetism, Wollan and Shull

Rietveld Analysis refines atom positions using powder data

Moderate resolution and rate requirements

Very large detector area required – Cost is important!

Oxide Ion Conductor – $\text{Sr}_8\text{Fe}_8\text{O}_{23}$



Correct Crystal Structure




Fig. 5

Single Crystal Diffuse Scattering

Bragg peaks give crystal structure

Light atom positions

→ especially *Hydrogen*

Magnetic structure

Diffuse scattering gives information on a larger length scale, such as:

Magnetic short range order



Hydrogen distribution in metals

Contrast between metals in alloys

Example: Cubic Stabilized Zirconia

Bragg peaks are several orders of magnitude larger than the diffuse scattering, but there is little separation in reciprocal space

→ need *resolution* and *dynamic range*



Proffen et al. (1996). *Acta Cryst.*, **B52**, 66-71.

Detector Needs for Diffraction



“Large-Pixel” (powder, engineering, disordered)

Very large detector area – low cost required

High efficiency requirement

20% at 50 eV (disordered)

50% at 1 eV (powder, engineering)

“Small-Pixel” (single crystal)

High resolution → 1 mm x 1 mm or better

Large area, low cost needed

Low γ sensitivity

Minimal magnetic effects for magnetic or polarized neutron work

Engineering imaging device – 100 micron resolution

Position Sensitive Monitors for Incident Beam

Large Scale Structure Diffraction



Small Angle Neutron Scattering (SANS)

Small Angles \Leftrightarrow Low Momentum Transfer \Leftrightarrow Large Structures

Metallurgy, Geology, **Polymers !**, Biology

Nanoscale structures

Requires high resolution 2D detector

Low noise, large dynamic range

Neutron reflectometry \Leftrightarrow 1D probe of surfaces and interfaces

“Sees” scattering length density differences between layers

Magnetic structure of films

Manufactured multilayer structures


Polymers, **complex fluids**, biological membranes

Interfaces (inside larger samples)

Detectors similar to SANS

\Leftrightarrow Total reflection is important \Leftrightarrow **Very High Rates !**


SANS and Reflectometer Instruments




KWS2 SANS at Jülich



Riso SANS at PSI



NG1 Reflectometer at NIST



Liquids Reflectometer
at SNS



Neutron Reflectometry




Total reflection $> 10^6$ n/s/1mm² pixel (10^8 n/s)

Dynamic range 1 to 10^{-10}

High Rate, High Resolution,
Low Noise, Low Gamma Sensitivity



Melittin in alkanethiol/phospholipid
hybrid bilayer membranes



Krueger *et al.*, *Langmuir* **17**, 511 (2001)
Majkrzak *et al.*, *Biophys. J.* **79**, 3330 (2000)

Detector needs for SANS & Reflectometry



In many cases, detectors require attenuators for low angle positions, not only to prevent saturation but to avoid detector damage.

In many research areas, (melting, self-assembly of membranes, surface/interface formation) there is a desire to observe time-dependent processes on the shortest possible time scale – maybe $< 1\text{s}$!

→ *Not compatible with mechanical attenuators*

Improved spatial resolution = better instrument resolution

Better angle definition also helps resolution

Big samples, high fluxes – need low γ sensitivity

Spatially resolved incident beam monitor

Inelastic Scattering

The change in the velocity of the neutron reveals the momentum and energy of the excitation.

Phonons

Spin Waves

Quasielastic excitations

Molecular vibrations

$\sigma_{\text{inelastic}} \ll \sigma_{\text{elastic}}$

Small Signals

Large Samples

Crystal Analyzer Spectrometers:



Triple Axis, Backscattering

HYSPEC, SNS Backscattering





Direct Geometry Spectrometers:

Fermi Chopper, Multi-Chopper

ARCS, CNCS, SEQUOIA



Inelastic Instruments



Fermi Chopper Spectrometer (Direct Geometry)



Fermi Chopper Spectrometers access excitations from
~10 to ~500 meV over a large region in momentum space


Position Resolution **1 – 2 cm**

Time Resolution **1 μ s**

Large Samples,
Powder or Single Crystal

Count rates are modest,
need low noise,
low γ sensitivity

*But Bragg peaks can mean
high rates on small spots!*



D. A. Tennant *et al.* Phys. Rev. B 52, 13368 (1995)

Detector Needs for Inelastic Scattering



Present ^3He LPSD gas detectors do a good job of gamma rejection and have very low noise

Limitations:

Peak Counting Rates -- The inelastic signals have very low counting rates, but elastic single-crystal Bragg peaks can have count rates of $> 10^7$ n/s (more for triple axis), which can saturate detectors and obscure inelastic features

Time Resolution – Uncertainty in time of arrival and flight path distance affect energy resolution, especially with cold neutrons. In other words, a flat detector is desirable.

Spatial Resolution – Present inelastic instruments are designed for detectors with resolution around 2.5 cm. Better spatial resolution could lead to better momentum resolution, especially with small samples.