

Neutron Detector R&D Roadmap

Ron Cooper

SNS Detector Team Leader

May 29-30, 2003

Roadmap Components




- Science case
 - Presented by Lowell Crow
- Instrument requirements
 - 24 beam lines in high power target station
 - 14 approved instruments
- Detector deficiencies
 - Almost all instruments need better detectors
- Detector R&D
 - SNS driven at present
 - Prototype development
 - Education and outreach

Roadmap Components (2)



- Management plan
 - Maintain science focus
 - Communication and coordination
 - Organize pool of expertise, provide reviews and direction
 - Maintain whitepaper: http://www.sns.gov/documentation/Neutron_Detector_White_Paper_March_03.pdf

Instruments



Component 2: Instrument Requirements



Table 3. Instrument requirements

Instrument	Number of pixels	Pixel area (cm ²)	Maximum neutron energy (eV)	Neutron capture efficiency %	Gamma efficiency	Time resolution (μs)	Peak pixel count rate (n.s ⁻¹)	Detector count rate (n.s ⁻¹)	Data transfer rate (M b/s)
Powder Diffractometer	40,000	2.4	0.33	50	10 ⁻⁶	1	100	3.5 × 10 ⁶	28
Disordered Materials Diffractometer	150,000	0.25	50	20	10 ⁻⁶	1	300	4.2 × 10 ⁷	340
High-Pressure Diffractometer	100,000	0.02	0.5	50	10 ⁻⁷	1	1 × 10 ⁴	3.0 × 10 ⁵	2.4
Engineering Diffractometer	80,000	1.25	0.15	50	10 ⁻⁶	1	2 × 10 ⁵	2.4 × 10 ⁶	20
Single-Crystal Diffractometer	5 × 10 ⁶	0.01	0.35	50	10 ⁻⁶	10	2 × 10 ⁴	3.0 × 10 ⁵	2.4
SANS Diffractometer	40,000	0.25	0.08	50	10 ⁻⁷	10	1,500	2.0 × 10 ⁷	160
Liquids Reflectometer	40,000	0.01	0.02	50	10 ⁻⁷	10	1 × 10 ⁶	7.0 × 10 ⁷	560
Magnetism Reflectometer	40,000	0.01	0.03	50	10 ⁻⁷	10	1 × 10 ⁶	9.0 × 10 ⁷	720
Backscattering Spectrometer	4,500	1.3	0.01	50	10 ⁻⁶	1	1 × 10 ⁴	1.3 × 10 ⁵	1
ARC Spectrometer	70,000	2.5	1.0	50	10 ⁻⁷	1	1 × 10 ⁶ (Bragg)	5.0 × 10 ⁵	4
CNC Spectrometer	15,000	6.3	0.05	50	10 ⁻⁷	1	1 × 10 ⁶ (Bragg)	7.0 × 10 ⁶	56
HRC Spectrometer	70,000	2.5	1.0	50	10 ⁻⁷	1	1 × 10 ⁶ (Bragg)	4.0 × 10 ⁵	3.2

May 29-30, 2003

Critical Requirements



- Powder Diffractometer
 - Large area coverage $> 10 \text{ m}^2$ when completed
 - 6 mm x 40 mm pixels
 - 50% efficiency for 325 meV, 0.5 Å neutrons
- Disordered Materials Diffractometer
 - Large area coverage $> 5 \text{ m}^2$ when completed
 - 20% efficiency for 50 eV neutrons
 - Long term stability
- High Pressure Instrument
 - Low gamma sensitivity
 - Position resolution: 1.5mm x 1.5mm
 - Continuous coverage

Critical Requirements (2)



- Engineering Instrument
 - Large area coverage $> 10 \text{ m}^2$ when completed
 - Maximum pixel rate: 20 n in $100 \mu\text{s}$ from elastic peaks
 - $100 \mu\text{m}$ resolution residual stress detector
- Single Crystal Diffractometer
 - 1-mm resolution
 - High dynamic range with minimal pixel cross talk
 - Maximum rate: 2×10^4 N per second per peak, 2000 peaks
- SANS Instrument
 - Maximum rate: 2×10^7 N per second
 - Low gamma sensitivity
 - 50% efficiency for 80 meV neutrons

Critical Requirements (3)



- Liquids and Magnetism Reflectometers
 - 1-mm position resolution
 - Maximum rate: 10^6 N per second for 100 pixels, 10^8 total
 - High magnetic fields
 - Low gamma sensitivity
- Backscattering Spectrometer
 - Detectors in vacuum
 - Low gamma sensitivity
 - Stability
- CNCS Spectrometer
 - Large area coverage $> 20 \text{ m}^2$ when completed
 - Minimized TOF uncertainty
 - Recovery from Bragg peaks
 - Stability

Critical Requirements (4)



- ARCS and Sequoia Spectrometers
 - Large area coverage, $\approx 20 \text{ m}^2$
 - Detectors in vacuum
 - Recovery from Bragg peaks
 - 50% efficiency at 1 eV
 - Stability
- Fundamental Physics – part of experiment
- Triple Axis – He LPSDs
- Spin Echo – 30cm x 30cm at 3×10^7 N per second
- Protein Crystallography
- Chemical Spectroscopy
- Polarized Neutron Instrument

Component 3: Detector Deficiencies



Table 4. Detector deficiencies for SNS instruments

Instrument	Parameter	Desired	Current	Comment
Liquids & Magnetism Reflectometers	Pixel area (cm ²)	0.01	0.02	0.02 is state of the art for ³ He gas detectors
	Maximum instantaneous rate/pixel (counts/s)	1.3 × 10 ⁶	7 × 10 ⁴	Beam attenuator will be necessary
	Maximum total instantaneous rate (counts/s)	1.2 × 10 ⁸	1 × 10 ⁶	Beam attenuator will be necessary
	Maximum time average rate/pixel (counts/s)	6.2 × 10 ⁵	7 × 10 ⁴	Beam attenuator will be necessary
	Maximum total time average rate (counts/s)	5.9 × 10 ⁷	5 × 10 ⁵	Beam attenuator will be necessary
	Transmission monitor pixel area (cm ²)	0.04	10	Characterize angular dependence of inc. beam
Powder Diffractometer	Neutron efficiency at 1 eV (%)	50	30	60% reduction in data rate
	Detector cost (\$/m ²)	150K	250K	Wavelength shifting modules will cover more area for the same cost
	Transmission detector maximum time average data rate (counts/s)	3.4 × 10 ⁷	1 × 10 ⁶	Reduce uncertainty in beam normalization
Engineering Instrument	Spatial resolution (mm)	0.1	1.0	Needed for residual stress depth profile measurements
	Transmission detector maximum time average data rate (counts/s)	5 × 10 ⁷	1 × 10 ⁶	Beam attenuator will be necessary
Single-Crystal Diffractometer	Spatial resolution (mm)	1	3	Unit cells limited to 30Å or less
	Transparent scintillator brightness (photons/neutron)	30,000	10,000	Needed for 1-mm resolution detectors
	Dynamic range (peak counts/background counts)	1 × 10 ³	1 × 10 ²	Needed for diffuse scattering studies
Inelastic Chopper Spectrometers	Spatial resolution (mm)	10	25	Q resolution limited by detectors for small samples
	Time resolution (μs)	1	5	Needed for high-resolution energy measurements
	Maximum instantaneous rate per detector (counts/sec)	2 × 10 ⁷	7 × 10 ⁴	Detectors will saturate, and inelastic data will be lost
Disordered Materials Diffractometer	Detection efficiency for 50eV neutrons (%)	20	5	Needed to measure atomic connectivity and defect distributions
Extended-Q SANS	Maximum total time average rate (counts/sec)	5 × 10 ⁷	5 × 10 ⁵	Needed to study weakly scattering biological samples
	Maximum parallax error (mm)	5	20	Q resolution limited by detector parallax

Critical Deficiencies



Instrument	Rate	Resolution	Efficiency
Liquids & Magnetism Reflectometers	x100	x2 Spatial (x250 Monitor)	
Powder Diffractometer	(x35 Monitor)		x1.6
Engineering Instrument	(x50 Monitor)	(x10 Residual Stress Detector, 1-D)	
Single Crystal		x10 Spatial	x3
Inelastic Chopper Spectrometers	x300 (Bragg peak)	x5 Time x2.5 Spatial	
Disordered Materials Diffractometer			x5
EQ-SANS	x100	x4 (Parallax)	

Component 4: Detector R&D



- No Standard Detector
 - Size varies from 4 cm² to 20 m²
 - Efficiency requirements vary by 2 orders of magnitude
 - Maximum pixel rates vary by 4 orders of magnitude
 - Dominant requirement that drives detector selection can be:
 - Efficiency
 - Gamma sensitivity
 - Time resolution
 - Cost
 - Rate capability
 - Stability
 - Other?

- 2-D Gas Detectors needed for:
 - Low gamma sensitivity and high stability applications
 - New high rate applications lead to
 - Discrete pixels – 40,000 channels
 - Ionization mode ^3He
 - Multilayer GEM
 - Micromegas
 - Potential Applications
 - SANS Instruments
 - Spin echo Instruments
 - High pressure Instruments
 - Reflectometers

- Scintillator Detectors needed for:
 - Large area coverage with small pixels
 - High neutron detection efficiency
 - New bright, transparent scintillator
 - 1- μ s time resolution
 - New readout schemes
 - Wavelength shifting fiber
 - Standard fiber
 - Anger camera
 - Potential Applications
 - Powder diffractometers
 - Engineering Instruments
 - Single Crystal Diffractometers
 - Disordered Materials Diffractometers

- Proportional Counters needed for:
 - Low-cost large-area coverage
 - Stability
 - Low gamma sensitivity
 - New designs needed
 - Improved time resolution
 - Eliminate Bragg peak saturation
 - Potential applications
 - Backscattering spectrometers
 - Chopper spectrometers
 - Triple axis spectrometers

High Resolution Detector R&D




- High resolution detectors needed for:
 - Sub millimeter time dependent imaging
 - Small samples
 - Very high rate applications
 - Neutron beam characterization
 - Potential detector Research includes:
 - Semiconductor detectors
 - Microsphere detectors
 - Micromegas Detectors
 - Conversion foils
 - Potential Applications:
 - Residual stress measurements
 - Small samples on the High Pressure instrument
 - Transmission detectors

Component 5: National Plan for Neutron Detector Development



- Make sure neutron detector development meets the needs of the neutron users community
- Coordinate detector R&D within the U.S. and worldwide
- Prioritize the use of limited R&D resources
- Plan components:
 - Executive committee
 - Advisory panel
 - Neutron users (6)
 - Detector experts (6)
 - Roadmap
 - Whitepaper

Management Plan




Executive Committee



- Current membership:
 - John Cameron – IUCF
 - Pat Gallagher – NIST
 - Veljko Radeka – BNL
 - Christine Hoffmann – SNS
 - Ron Cooper – SNS
- Responsibilities:
 - Supervise the plan
 - Maintain communication
 - Administer workshops
 - Periodic reviews
 - Request support from the advisory panel

- Diverse groups



Roadmap Components



- Advisory panel
 - International membership
 - Pool of expertise to support
 - Funding agencies
 - Executive committee
 - Detector PI's
 - Review proposals
- Roadmap
 - Whitepaper
 - Reference for proposals
 - Living document
 - Expand to include steady-state sources
 - Update as needed

Summary



- The roadmap consists of:
 - Science case
 - Instrument requirements
 - Detector deficiencies
 - Detector R&D
 - Management plan
- Maintained as a whitepaper
 - Living document
 - Reference for proposals
- Goal is a coordinated effort
 - Science focus
 - Communication
 - Guidance

