

Neutron detectors development

Toward future spallation sources

Bruno Guerard

- MicroStrip gas chambers (MSGC)
- MultiWires gas chambers (MWPC)
- Position sensitive counter tubes

Developments supported by the ILL Millennium program

- large angle coverage banana MWPC
- The Multitube detector
- high resolution 2D detectors (scintillators and He3)


International scientific co-operation

The ILL is operated by its three founder countries:

France	31%
Germany	31%
United Kingdom	26%

In association with its partner countries :

Austria:	1.5%
Spain:	3%
Italy:	3%
Czech Republic:	0.5%
Russia:	2%
Switzerland:	3.5%



ILL/SDN (Service des Détecteurs de Neutrons)

Jean-Claude Buffet **Mechanical design**

Gilbert Viande **Mechanical design**

Michel Gamon **Electronician**

Fabrice Horst **Electronician**

Frédéric Millier **Electronician**

Jean-Francois Clergeau **Physicist**

Bruno Guerard **Physicist**

Giuliana Manzin **Physicist**

Patrick Van Esch **Physicist**

Fabrication outlines

Gas filling

outgasing, detection efficiency, position resolution

Electronics

individual readout versus charge division

Type of detector

MSGC, MWPC, PSCT, scintillators



Nuclear reactions for neutron detectors


- $n + {}^3\text{He} \rightarrow {}^3\text{H} + {}^1\text{H} + 0.764 \text{ MeV}$
- $n + {}^6\text{Li} \rightarrow {}^4\text{He} + {}^3\text{H} + 4.79 \text{ MeV}$
- $n + {}^{10}\text{B} \rightarrow {}^7\text{Li}^* + {}^4\text{He} \rightarrow {}^7\text{Li} + {}^4\text{He} + 0.48 \text{ MeV } \gamma + 2.3 \text{ MeV (93\%)}$
 $\rightarrow {}^7\text{Li} + {}^4\text{He} + 2.8 \text{ MeV (7\%)}$
- $n + {}^{155}\text{Gd} \rightarrow \text{Gd}^* \rightarrow \gamma\text{-ray spectrum} \rightarrow \text{conversion electron spectrum}$
- $n + {}^{157}\text{Gd} \rightarrow \text{Gd}^* \rightarrow \gamma\text{-ray spectrum} \rightarrow \text{conversion electron spectrum}$
- $n + {}^{235}\text{U} \rightarrow \text{fission fragments} + \sim 160 \text{ MeV}$
- $n + {}^{239}\text{Pu} \rightarrow \text{fission fragments} + \sim 160 \text{ MeV}$

Advantages of ${}^3\text{He}$ detectors

- counting dynamics
- Time resolution
- low sensitivity to γ

3He Pulse Height spectrum







Increasing the CF_4 pressure reduces the difference between the centre of gravity of the primary charge cloud and the point where the neutron interacted.

2 bars CF_4
1.4 mm FWHM

Detection efficiency and scattering



Scattering & absorption measurement performed with a position sensitive counter tube and a sample of Aluminium (5083) 5 mm thick
 measured attenuation : 5.92% (abs. + scatt.)
 Calculated absorption : 1.72 %
 ⇒ scattering = 4.2%





Neutron conversion efficiency

$$\varepsilon = 1 - e^{-gap \cdot \sigma \cdot n}$$
 gap = 3 cm; n = density of gas
 s = absorption cross section
 = 5333 barns * $\lambda/1.8$

Gas purity

Sealed detectors --> require very pure gas

- Selection of the mechanical components of the detector
- Double O-Ring, or metallic gasket
- heating at 150°C under control of a mass spectrometer.




MSGC detectors

Advantages :

- Spatial uniformity

Drawbacks :

- One supplier only (IMT)
- Limited size



D20


Angular coverage : $160^\circ \times 5,8^\circ$

Angular resolution : $0,1^\circ$

working since July 2000


D4C



Modularity & high efficiency





Gabriel Ouello
www.jll.fr/YellowBook/D4C/

Aperture : 145° (2 scans) \times 10 cm height
 Position resolution : 2.5 mm
 Absorption gap : 30 mm
 Gas pressure : 15 bars ^3He + 0.3 bar CF_4
 Detection efficiency : 90% (@ 0.7 angstroms)
 Counting rate : 50 kHz/cell








*Response with a vanadium sample
Maximum deviation from average :
5% before calibration
1% after calibration*

*Nickel sample diffraction curve
used to calibrate the position of
the 9 MSGCs*




Diaphragm

MSGC

80 mm × 80 mm × 1 mm

5 detectors working at ILL
the first without interruption since 96


86 mm



Rear side



Charge division readout

Individual Readout 64x64
Position resolution : 3 mm
Useful area: 192 mm x 192 mm
94% efficiency at 2.4 Angstroms
Counting rate : 500 kHz (with
base line restorer amplifiers)



*TiO₂ powder sample
@ 2.5 Angstroms*





Fabrication features

Substrate : S8900 240 mm x 240 mm x 1.5 mm.

Entrance window : 5 mm aluminium

Gas Mixture : 5 bars ^3He + 2 bars CF_4

working since August 2001

Virtual cathode layout

Higher maximum amplification factor --> more robust

Anode and cathode signals are equal (and of opposite sign)

Double connection of the anodes

- > parallel charge division readout


- > Easy control of anode resistance value

But counting rate limitation


Anodes



Cathodes



X-Rays



^{55}Fe source pulse height spectrum in Argon:CH₄ (90:10)
(source inside the detector vessel)


$$G = 10^5$$

XY individual readout : Dead time & counting limits

1 μ s/strip (time development)


100 ns global (XY coincidence)

3 KHz/mm of anode (glass charging)





XY charge division along each anode strip


- Improved position resolution in one direction
- No XY coincidence dead time




Charge division along the anodes : position linearity



Position resolution



1.3 mm FWHM = limit of the
stopping gas (2 bars CF_4)




Future work : improve the position resolution and the counting rate


- divide by 5 the substrate thickness (requires a gluing of the glass on a metallic reinforcement) --> 0.3 mm
- multiply by 5 the number of anodes --> 0.6 mm pitch
- add floating electrodes between anodes to help discharging the glass surface

MWPC detectors

Bidim26 MWPC



working since February 2002



Bidim26 MWPC





PSCT


(Position sensitive counter tubes)



Actual D22 detector

2 orthogonal frames of sensing wires with a pitch of 7.5 mm.

A coincidence time of 200 ns between the 2 frames limits the counting performance of this type of detector to about 500 kHz (200 kHz with all contributions)



SANS_2MHz

By Measuring the signals on both ends of a resistive wire centred in a tube, the position of each neutron can be determined along the tube without affecting the measurement on the other tubes. This principle has been adopted to improve the counting rate on the future D22 detector.

Small diameter PSCT

8 mm diameter, 1 m long, and 15 bars of He3.




Electronics

Low noise amplifier, pulse shaping, base line restorer circuit, followed by a charge division processing card which converts the two analogue signals into the position.


Position resolution = 6 mm (3 mm at very high gain)

Counting rate


From the measured rate versus true rate relationship, the dead time of one PSCT associated with its electronics is estimated at 1.4 μ sec. If we define a limit of 10% counting loss, the maximum counting rate of one PSCT is then 70 kHz. For an uniform irradiation, this value corresponds to 9 MHz for the future assembling of 128 PSCT, a factor of 50 better than the actual detector.



efficiency along detector



Specification of spatial resolution comfortably met :
6 mm (9 mm at 100 KHz/tube).



A pilot production of 16 PSCT will be tested with the final version of the electronics and of the mechanics

In order to avoid High voltage sparks, the mechanical support must guaranty the anode wire centring with an accuracy of 0.1 mm along the tube length.

MAPS - 16 m² with 15 x 25 mm resolution




The MAPS spectrometer at ISIS is optimised to measure high energy magnetic excitations ($> \sim 50$ meV) in single crystals with energy resolution 1-5% depending on choice of monochromating chopper. The principal innovation in MAPS is the use of position sensitive detectors that provide close to continuous coverage over a large solid angle detector array in the forward direction.

Applications

- Magnetic excitations in single crystals
- Structural excitations in single crystals
- High resolution studies of non-dispersive excitations

Detectors

- 6 m from sample position
- approx 16m² array of 147,456 pixel elements
- 574 Position sensitive ³He tubes (10 bar partial pressure) 2.5 cm diameter, resolution 15mm along their lengths







The Multitube detector

**x tubes are brazed on
a common gas vessel
on both ends.**


Fabrication of the MT32 prototype









Position resolution




The MT4 prototype

4 tubes of 1 m long and 8 mm diameter.




**measured detection
parameters = those obtained
with PSCT.**

Number of counts




Anode voltage

Position resolution FWHM (mm)




Spatial linearity

Position resolution




The Multitube as an alternative for the detector of the IN5 TOF spectrometer


10 curved modules, 32 tubes each, 2 m high.




*A very large angular
coverage banana MWPC*





Assembling of 40 segments on one cradle



View of one segment






Cathode strips 1 mm wide

Interlinked with a resistive line of ~ 10 kohms

Substrat : Sodalime glass 410 mm x 40 mm x 4.8 mm

Chromium metallisation



Specifications

Radius of curvature (sample to conversion gap) : 73 cm

Active area : 40 segments of (4 cm x 40 cm)

Angular coverage : $30^\circ \times 120^\circ$


Spatial resolution : 2.5 mm in X (0.19°), 2 mm in Y (0.15°)

Gas / Conversion Gap : 4.5 bars He^3 / 3 cm

Detection efficiency @ 1 A : 60%

Total readout channels : 640 in X, 40 in Y

Counting rate : 1 MHz global and 50 kHz local



Calculation results

Maximum Strain values


vessel : 1.21 mm

Window : 0.68 mm

Maximum Stress value

calculated on the window

37 Mpa




(Strain scale multiplied by 200)

A very low efficiency beam monitor

Convertor	N ₂	³ He
Nuclear reaction	$^{14}\text{N} + \text{n} \rightarrow \text{p} + ^{14}\text{C} + 627 \text{ keV}$	$^3\text{He} + \text{n} \rightarrow \text{p} + ^3\text{H} + 764 \text{ keV}$
Capture cross section @1.8 Angstroms (mbarns)	1.8	5350
Gas pressure to obtain 10 ⁻⁵ detection efficiency in 1 cm conversion gap	100 mbars	0.07mb

Pulse Height spectrum

gas mixture : 12 mb N_2 + 400 mb CF_4





Scintillators

ZnS:Ag/Li6F scintillator

Thickness : 0.42 mm

Decay to 10% : **80 μ s**

Emission spectrum peaked at 450 nm
~ 160.000 photons/ neutron




Li Glass scintillator (GS20)



Thickness : 0.5 or 1 mm

decay to 10% : **60 ns**

Emission spectrum peaked at 395 nm
~ 7.000 photons/ neutron



Scintillator : ZnS:Ag/Li6F
0.5 thick square WLS fibers Y11
4x4 multi-anode PM



Position sensitive PM : Scintillator directly in contact with the PM

- + : cost and gamma discrimination
- : counting rate






Image measured with a Cadmium mask in front of the scintillator-PM detector.

Position resolution : 0.6 mm (FWHM)



Scintillator	Neutron detection ϵ_n	^{60}Co Gamma sensitivity ϵ_γ
$^6\text{LiF/ZnS:Ag}$	52 %	$5 \cdot 10^{-8}$
GS20 (0.5 mm)	60 %	$2 \cdot 10^{-7}$
GS20 (1 mm)	78 %	$1 \cdot 10^{-5}$



GEM - fibre coupled zinc sulphide scintillator detectors

GEM is designed to be a high intensity, high resolution neutron diffractometer for structural studies of disordered materials and crystalline powders.

The detectors will cover a very wide range in scattering angle from below 5° to 170° .

When GEM is completed ...

~3.5 steradians



Area 10 m²

6550 detector elements in 75 modules.

660,000 individual optic fibres, whose total length will be about 350 kilometres.





X resolution 5 mm
 Y resolution 200 mm
 detector efficiency (^{60}Co) 10^{-7}
 intrinsic background 12 c/element/hr
 pulse pair resolution $2.5\ \mu\text{s}$
 stability $< 0.1\ \%/^{\circ}\text{C}$




120 elements
 2C_n coding to 16 PMT
 coincident pairs of PMT outputs determine
 origin of a single neutron interaction



Active gas Scintillators (GEM + CCD camera)



Visible light produced by
interactions of neutrons





- ✦ The reliability of MSGC detectors is now established. They do not provide a significant improvement as compared to MWPC.
- ✦ For continuous sensitive area larger than 20 cm x 20 cm, good position resolution is obtained with MWPC.
- ✦ PS Counter tubes, and similar detectors, will be extensively used in SNS for several large area instruments with moderate position resolution. The Multitube provides an effective way to decrease the fabrication cost.
- ✦ Scintillator detectors with counting or integrating devices are an alternative solution to gas counters