

In-Situ Formation and Alignment of Biological Supramolecular Assemblies in Microchannel Devices

Youli Li

Materials Research Laboratory
University of California, Santa Barbara

In Collaboration with (at UCSB)

The Safinya Group: Uri Raviv, Dan Needleman, Linda Hirst,
Cyrus Safinya, ...


The Wilson Group: Hurbert Miller, Les Wilson

The MacDonald Group: Emily Parker, Noel MacDonald


Funding Support: NIH GM-59288, NSF DMR- 0203755, CTS-0404444, DOE W-7405-ENG-36 (BES)

Alignment of self-assembled structures in microchannels

- Design a system to **geometrically confine** filamentous biomolecular self-assemblies.
- we must take into account the persistence length of the molecule/structure.



- Create a large **array** of confined molecules, to produce a highly aligned protein sample.




- Use **x-ray diffraction** to study small amounts of aligned protein samples

Other alignment methods – shearing, magnetic fields
larger amounts of protein


Microchannel system Design

- The micro channels are designed to be easily filled by a user by capillary flow
- The channels are coated with a PEG 'comb' to resist protein absorption.
- The wells are used to fill the channels with different protein solutions.
- Bundles self-assemble in the channels



Microfabrication of the microchannel

Silicon




Standard
BOSCH etch
process

High aspect
ratio features


Channel
structure
etched into
Silicon wafer

SiO₂ thermally
grown on the
surface

PEG-Silane
absorbed to
the surface



Titanium




Chlorine based
TIDE etch.
High aspect
ratio.
Smooth walls
robust material
Thin and
strong

Channel
structure
etched into
Titanium foil

TiO₂
Sputtered on
the surface


PLL-PEG
absorbed to
the surface



Filling the Microchannels with actin/ α -actinin bundles


- G-Actin is polymerized *in-situ* in the presence of the linker molecule - Like behavior in the cell.
- Bundles grow along the channel, aligned by confinement.
- The wider the channel, the less alignment is observed.

α -actinin conc. \longrightarrow




Bundle growth in microchannels


Confocal fluorescence microscope images



No channels




10 μ m Ti channels




20 μ m Ti channels

Microchannel Device X-ray Scattering Results



Intensity (a.u.)



•Samples prepared in the channels system give **highly aligned** results with **very little material** compared to a typical capillary sample.

•For our **actin/ α -actinin system** we observe the expected peaks for a quasi-square structure.

N. Bouxsein, L.S. **Hirst**, Z. Abu Samah, Y. Li and C.R. Safinya *APL* **85**, 5777 (2004)

L.S. Hirst et al. *Langmuir* **21**, 3910 (2005)

X-ray & Neutron Requirements

- Micro- to nano-sized beams
- In-situ multi-probes at multiple lengthscales
- Parallel automated sample handling and data acquisition
- Issues with signal and radiation damage