



CHESS & LEPP

Impact of a Future Energy Recovery LINAC X-ray Source on Nanoscale Science

Joel D. Brock*

*School of Applied and Engineering Physics
Cornell University, Ithaca, New York
jdb20@cornell.edu*

* on behalf of the ERL team of I.V. Bazarov, S.Belomestnykh, D. Bilderback, J.D. Brock, E. Fontes, K.D. Finkelstein, S.M. Gruner, G. Hoffstaetter, A.Kazimirov, M. Liepe, Y.Lin, H.S. Padamsee, D.Sagan, V.Shemelin, Q. Shen, K.Smolenski, C. Sinclair, R. Talman, M. Tigner, V.Veshcherevich, Cornell University; and L. Merminga & G. Krafft of J. Lab



Energy Recovery LINAC




CHES & LEPP

Injector: high-brilliance electron bunches generated by a laser-driven photocathode are accelerated to $\sim 10\text{MeV}$

Main Linac: superconducting cells accelerate electron bunches to 5 GeV, and recover energy from returning bunches (180° out of phase)

Transport loop: produces high-brilliance x-ray beams through undulators, and re-injects electrons into main linac for energy recovery





- Accelerating bunch
- Returning bunch



Electron Beam Size/Shape



CHESS & LEPP





Advanced Photon Source compared with Energy Recovery Linac



CHES & LEPP

Parameter	APS 3rd generation storage ring	Energy recovery linac	Gain factor
Electron source size in microns rms	239(h) x 15(v)	2(h) x 2(v)	1/900 in area
Micro x-ray beamsize	100 nm to 1 micron (planned 30 nm)	1 nm	100 to 1000
Coherent flux x-rays/sec/0.1% bw	3×10^{11}	9×10^{15}	3,000
Pulse duration (rms)	32 ps	<100 fs	more than 320 times shorter

Conclusion: the ERL machine will be transformational!




ERL Will Provide Unprecedented Nanobeams



CHESSE & LEPP

Storage ring nanobeam flux limited by source size, shape, and divergence.



- Intense 1-10 nm probe size (rms), 1-10 keV beam allows study of nanostructures and molecules
- Quantitative atomic-scale structure, strain, orientation imaging
 - Increase fluorescent trace element sensitivity from present 10^{-19} g to single atom (10^{-24} g)
- Sensitive to chemical state via XAFS at ultra-low concentrations
- Ability to penetrate thick layers, nasty gas environments, etc. (as opposed to EM)


ERL source with electron beam size of 2 microns rms for 1 m long undulator and 0.5 m beta function demagnify by 2000x to make 1 nm beam size, etc.



Microbeam X-ray Science




CHESS & LEPP



Zinc distribution in plant leaf by
SR x-ray
fluorescence

few cm scale object (CHESS
data)


Centimeter
scale



Hot-rolled Aluminum

SR x-ray diffraction. Map grain
orientation and stress in real
samples of 10^4 cubic microns at 1
micron resolution (APS data)

Micron scale



Two impurity atoms (yellow dots)
in silicon crystal

TEM with 200 keV electrons
can see individual atoms on
samples a few atoms thick
(Voyles, Lucent Technologies)

Nanometer
scale

Atomic scale





Future of X-ray Nano-Characterization



CHESS & LEPP

Techniques

- Diffraction, Scattering, ...
- Spectroscopy: EXAFS, XANES, ...
- Holography
- Tomography
- Imaging: Phase contrast, Coherent, and Fluorescent (single atom sensitivity, $\sim 10^{-24}$ gm)

50 nm
nanocluster

1–10 nm Probe Size
(diffraction limited)

