

## Science Opportunities on Early Spallation Neutron Source Instruments

K. W. HERWIG

Instrument Scientist

and

Deputy Senior Team Leader for SNS Instrument Systems

**October 8, 2004** 



#### **The Spallation Neutron Source**






- The SNS will begin operation in 2006 in Oak Ridge, Tennessee
- At 1.4 MW it will be ~8x ISIS, the world's leading pulsed spallation source
- The peak thermal neutron flux will be ~20-100x ILL
- Initial facility will have 1 target station with 24 instruments
- An upgrade to 3-4 MW and a second target station is planned
- SNS will be a short drive from HFIR, a reactor source with a flux comparable to the ILL



#### **Instrument Suite**





# **Backscattering Spectrometer Installation**










# **Reflectometer Components**








# **POWGEN3** Components






#### **Instrument Layout**

SNS SING UNIV NUC-PH Other






OAK RIDGE NATIONAL LABORATORY U. S. DEPARTMENT OF ENERGY



## Timeline for scientific productivity








#### Phase Space probed by Inelastic Instruments









## Phase Space probed by Elastic Instruments










#### **Characteristics of Friendly Users**



- PATIENT Intermittent performance of Instrument/Target/Accelerator
- ENGAGED Interested in participating in debugging instruments/software



- MOTIVATED Getting the data may take more effort than usual
- LOW MAINTENANCE SAMPLES may need to stay in beam for extended periods of time, limited sample environment equipment in 2006

#### **EXCITING SCIENCE**



# Unique Features at Low Power/Duty Cycle (pre-January 2008)



- World class accelerator performance in Early 2008
- Instrument designs typically multiplied source power gains by factors of 2 to 10
  - Capable of world class science prior to 2008
- Early Capabilities
  - Backscattering Spectrometer
    - Wide dynamic range enables study of dynamics over energy transfer range of 2 µeV to ~18 meV (~ 4 decades) – glass/amorphous materials dynamics
  - Reflectometers
    - Flexible Data Acquisition System allows "pump-probe" techniques – application of external fields and measured response on ~ 1 minute time resolution
    - Time-of-flight coupled with area detectors facilitates offspecular reflectivity – in-plane structures, magnetic domain formation, nano-patterned films



# Unique Features at Low Power/Duty Cycle (pre-January 2008) #2




- Early Capabilities
  - ARCS (Fermi Chopper inelastic spectrometer)
    - Wide detector angular range/position sensitive elastic diffuse scattering
    - Density of states vs. T and P (large powder samples)
  - CNCS (Cold Neutron Chopper Spectrometer inelastic)
    - Hydrogen diffusion in membrane materials
    - Wide angular range/position sensitive detector elastic diffuse scattering
  - POWGEN3 (general purpose diffractometer)
    - High Resolution (∆d/d ~ 0.15%) ab initio structure determination
    - Fastest U.S. neutron powder diffractometer
    - Large samples parametric studies, T, P, pO<sub>2</sub>



#### 16 Instruments are formally approved







#### CONCLUSION



- Web sites for the Spallation Neutron Source
  - www.sns.gov project public site
    - Follow links to user information and to instrument data sheets www.sns.gov/users/instrument\_data\_sheets.pdf
    - Follow links to ORNL Neutron Scattering Points of Contact (HFIR and SNS) personnel http://www.sns.gov/users/points\_of\_contact.pdf
- Instrument Commissioning
  - 3 in 2006
  - 5 more in 2007
  - 2 more in 2008
- Late 2007 Early 2008: at ISIS MWatt \* Hours level
- Opportunity for Friendly Users as instruments go through commissioning
- User Program Experiments 2nd Half 2007

