

***In -situ*, High Temperature Study of Phase Transformations in Ceramics**

Professor W. M. Kriven

J. Schneider, L. F. Siah, W. Yoon, K. Jurkschat, P. Sarin

**Department of Materials Science and Engineering
University of Illinois at Urbana-Champaign
Urbana, IL, USA**



Supported by NSF- DMR and AFOSR




Goal of this work

- To explore the feasibility of using the thermal-image technique to study reversible, displacive-type phase transformations, specifically in air, at elevated temperatures
- To measure their lattice parameters as a function of temperature and determine their axial thermal expansion coefficients, transformation temperatures and unit cell volume/shape changes

Halogen lamp



Osram Xenophot
HLX 64635
15V 150W




Radiation of the halogen lamp


A Compact Thermal-image Furnace

Heat source: halogen infrared reflector lamp (OSRAM Xenophot HLX64635)


4 lamps (50mm diameter, 150W) arranged to have a common focus



Intersection of filament image → "hot-spot" at the common focus



High Temperature Furnace




Water cooling

Goniometer for
specimen rotation


Thermocouple

Setup at UNICAT beamline 33BM at APS;
Argonne National Laboratory


Layout




Temperature Profile of Quadrupole (at 2/3 power)



Temperature measurement




Temperature measurement
Hot spot with Al_2O_3 tube
(OD : 5 mm, ID : 4 mm)



Temperature measurement with
change of lamp and tubes

In-situ HTPXRD: RNbO_4 (R = Dy, Y)

YNbO₄-Pt 1700C NSLS Sep 2001




Obvious TDS indicating attainment of very high temperatures *in-situ* in air in this work

Preparation of Capillary Sample

- Powder preparation
 - PVA method
 - ethylene glycol method
- Capillary fabrication
 - extrusion
 - “controlled” sintering
- “Thermometer” coating
 - platinum ink
 - oxide layer

(CeO_2 mp: 2300°C , MgO mp: 2800°C)



A representative capillary sample
(diameter $\sim 0.35 - 0.40$ mm)

Temperature Calibration


Up to 1700 °C

- Thermocouple Pt 30% Rh - Pt 6% Rh in close proximity to specimen
- Accuracy ± 10 °C at 1500 °C

Above 1700 °C

- Lattice constant of coating material
- Furnace power calibration with known melting points and phase transformations
- Accuracy ± 30 °C at 2000 °C

Reflection Geometry for *in-situ* high temperature phase transformation studies at UNICAT 33BM Beam line at APS



Advantages

- In situ studies of lattice parameter development and phase transformations possible up to 2000 °C
- In air studies
- Quick temperature ramping
-
- Easy sample preparation

Problems

Temperature calibration

- Insufficient known calibration materials
- Thermal expansion mismatch between coating and specimen
- Reaction between calibration material and specimen under study





Specimen stability

- Structural (specimen disintegrating, glassy phase)
- Chemical stability

X-ray absorption by heavy elements in Debye-Scherrer geometry





Thermal expansion of sample holder in Bragg-Brentano geometry

In-situ HTPXRD in Air: RNbO_4 ($\text{R} = \text{Dy}, \text{Y}$)



DyNbO_4 : monoclinic \rightarrow tetragonal transformation


In-situ HTPXRD in Air: RNbO_4 ($\text{R} = \text{Dy}, \text{Y}$)





DyNbO₄: monoclinic → tetragonal transformation

In-situ HTPXRD in Air: RNbO_4 ($R = \text{Dy}, \text{Y}$)

YNbO_4 : monoclinic-to-tetragonal transformation



Temperature evolution of
a, b, c, γ and unit cell volumes



Conclusion for LnNbO_4

- the monoclinic-to-tetragonal transformation
 - in YNbO_4 is 2nd order at ~830 °C
 - the room temperature spontaneous strain of YNbO_4 is 6.35%

<u>Compound</u>	<u>Crystal Symmetries</u>	<u>Transformation Temperature (To on cooling)</u>	<u>Volume Change (ΔV)</u>	<u>Unit Cell Shape Change($^{\circ}$)</u>
ZrO ₂	Tetragonal → monoclinic	950	(+)4.9% (R.T.)	9
Ln ₂ O ₃ (type)	monoclinic → cubic	600—2200	(+)10%	10
Ca ₂ SiO ₄ (K ₂ SO ₄ -type)	monoclinic → orthorhombic	490	(+)12%	4.6
Sr ₂ SiO ₄ (K ₂ SO ₄ -type)	orthorhombic → monoclinic	90	0.2%	2
NiS	rhombohedral → hexagonal	379	(+)4%	—
2Tb ₂ O ₃ .Al ₂ O ₃ (type)	orthorhombic → monoclinic	1070	(+)0.67%	18.83
PbTiO ₃	cubic → tetragonal	445	(+)1%	0
KNbO ₃	tetragonal → orthorhombic	225	~0%	0
LuBO ₃	hexagonal → rhombohedral	1310	(+)8%	—
MgSiO ₃ (CaSiO ₃ -type) (FeSiO ₃ -type)	orthorhombic → monoclinic	865	(-)5.5%	18.3
YNbO ₄ (LnNbO ₄ -type)	tetragonal → monoclinic	900	(-) 1.8%	4.53
LnBO ₃ (type)	hexagonal → hexagonal	550—800	(-)8.2%	—

Table 1. Examples of First Order Displacive Transformations in Ceramics


<u>Compound</u>	<u>Crystal Symmetries</u>	<u>Transformation Temperature (To on cooling)</u>	<u>Volume Change (ΔV)</u>	<u>Unit Cell Shape Change($^{\circ}$)</u>
Cristobalite (SiO ₂)	cubic → tetragonal	265	(-) 2.8%	0
Hexacelcian (BaAl ₂ Si ₂ O ₈)	hexagonal → orthorhombic	300	(-) 0.43%	0
Leucite (KAlSi ₂ O ₆)	cubic → tetragonal	620	~0	0
Zircon (ZrSiO ₄)	monoclinic → tetragonal	827	?	?
Di-lanthanide aluminates (Ln ₄ Al ₂ O ₉)	monoclinic → monoclinic	1400	(+) 0.5%	?
Di-lanthanide titanates (Ln ₂ TiO ₅)	hexagonal → orthorhombic	1712	?	0
Barium orthotitanate (Ba ₂ TiO ₄)	?	?	?	?
Cerium pyrosilicate (CeSiO ₄)	?	?	?	?
Aluminum titanate (Al ₂ TiO ₅)	?	?	?	?
Lithium phosphate (Li ₂ PO ₄)	?	340	?	?
Lanthanide (eg.Gd) vanadates (LnVO ₄)	monoclinic → tetragonal	825	?	?

Table 2. Other Examples of Phase Transformations in Ceramics

Table 3. RARE-EARTH OXIDES : PHASE TRANSFORMATIONS

Compound	La	Ce	Pr	Nd	Pm	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb
LnNbO ₄ scheelite	Monoclinic (low T) → Tetragonal (high T)													
Ln ₂ TiO ₅	Orthorhombic						Orthorhombic (low T) → Hexagonal → Cubic					Cubic		
LnAlO ₃ Perovskite	Rhombohedral (low T) → Cubic			Rhombohedral (low T) → orthorhombic → Cubic		Rhombohedral (low T) → orthorhombic		Orthorhombic						
LnTaO ₄	Unknown structure		Monoclinic (low T) → Tetragonal									Monoclinic		
LnVO ₄	Monoclinic (?)		Tetragonal YPO ₄ structure											
LnAsO ₄	Monoclinic CePO ₄ structure			Tetragonal YPO ₄ structure Tetragonal scheelite (CaWO ₄) structure (high pressure)										
LnPO ₄	Monoclinic CePO ₄ structure									Tetragonal YPO ₄ structure				
Ln ₂ Ti ₂ O ₇	Monoclinic			Cubic (pyrochlore)										
Ln ₃ NbO ₇	Tetragonal (?)			?			Hexagonal				?			


Toughening of Ceramics by Transformation Weakening of Interphases




Toughening Mechanism

- Schematic diagram illustrates “transformation weakening of ceramic interphases” leading to overall toughening of a ceramic matrix composite.
- In thermally induced transformations, all interphases are pre-transformed before the approach of a crack, with some consequent loss of overall strength of the material.
- In the ideal, shear-stress induced case, an oncoming crack induces a transformation in its immediate environment, with strength only minimally reduced throughout the bulk.
- Maximum toughening is achieved, since the propagating crack needs to do work to overcome the nucleation barrier and cause transformation, and onset of the other synergistic toughening mechanisms (e.g., crack formation) occurs.


Work of Fracture of TW Weakened Composite



Preferential Crack Deflection along TW Interphases




Toughening by TW of Interphases



Long term vision for study of phase transformations in structural ceramics and composites

Displacive phase transformations can be accompanied by volume changes and/or by shape changes

Large volume changes (ΔV)



- “Transformation toughening” (e.g. zirconia ceramics)
- “Transformation weakening” of interphases, giving overall toughening of fiber-reinforced composites or fibrous monoliths (e.g. enstatite, cristobalite)



Structural ceramics for airplane and
Ground-based turbine engines

Unit cell shape changes ($\Delta\beta$)



- Ferroelastic phase transformations
- Ferroelasticity
- Toughening by “twin flipping”



Large force actuation
(e.g., for MEMS devices)
Shape memory ceramics